



KEY NOTES

1. SEE ELEVATION OF A1107, UTILITY RAIL SYSTEM, FOR UTILITIES.
2. PROVIDE A1165, IV TRACKS, FOR NAVY PROJECTS ONLY.
3. PROVIDE AT LEAST ONE HANDWASHING FACILITY OR ONE FOR EVERY THIRD PREP/INDUCTION MODULE. P3200, ELECTRONIC LAVATORY, MAY BE USED INSTEAD OF P3100, LAVATORY. FOR THE P3200 PROVIDE EITHER BATTERY POWER OR HARD WIRED CONNECTION.
4. MOUNT A5145, HOOK, AT 1372 [4'-6"] AFF.
5. M7845, BEDSIDE PHYSIOLOGICAL MONITOR, MAY BE CONNECTED TO A M7850, CENTRAL PHYSIOLOGICAL MONITOR.

OPERATING ROOM PATIENT PREP./INDUCTION (ORPP1)

SURGERY (INPATIENT & AMBULATORY)

FLOOR/EQUIPMENT PLAN 37.25 NSM (401 NSF)

DATE: 01/28/2002

0 1000 2000



1:50 SCALE (MILLIMETERS)