



NOTES

1. ALL RECEPTACLES IN HEADWALL SHALL BE HOSPITAL GRADE.

KEY NOTES

1. PROVIDE A1165, IV TRACKS FOR NAVY PROJECTS ONLY.
2. VACUUM SLIDE.
3. SWITCH CONTROLS CEILING MOUNTED FLUORESCENT LIGHT FIXTURES.
4. TONE/VISUAL OR AUDIO/VISUAL AS REQUIRED BY USING MILITARY DEPARTMENT.
5. CODE BLUE AS REQUIRED BY USING MILITARY DEPARTMENT.
6. M7850, CENTRAL PHYSIOLOGICAL MONITOR, AND M7845, BEDSIDE PHYSIOLOGICAL MONITOR, MAY BE ADDED. VERIFY WITH USER.
7. P8705, DIALYSIS BOX, WITH COLD WATER AND DRAIN.
8. TWO LOW VOLTAGE MOMENTARY CONTACT TOGGLE SWITCHES TO CONTROL UP AND DOWN LIGHT SECTIONS OF FIXTURE M7435.



KEY NOTES

1. P8705, DIALYSIS BOX, WITH COLD WATER AND DRAIN.
2. PROVIDE A1165, IV TRACKS, FOR NAVY PROJECTS ONLY.
3. P3200, ELECTRONIC LAVATORY, MAY BE USED INSTEAD OF P3100, LAVATORY. FOR THE P3200 PROVIDE EITHER BATTERY POWER OR HARD WIRED CONNECTION. PROVIDE ONE P3100 AND ACCESSORIES FOR EVERY 2 DIALYSIS STATIONS.

General Notes:

1. For typical equipment mounting heights see "TYPICAL EQUIPMENT MOUNTING HEIGHTS" drawing.
2. MIL-STD-1691 data can only be accessed from the Equipment Plan.

RENAL DIALYSIS STATION, CHAIR (RDC01)

SPECIALTY MEDICAL CLINICS

ELEVATIONS

13.40 NSM (144 NSF)

DATE: 01/28/2002

0 1000 2000

1:50 SCALE (MILLIMETERS)